

Codonics DVB5 Hard Cassette Media Packaging & Storage

Technical Brief



Summary

This Technical Brief covers the date format, "Best By" date, storage and storage conditions for hard cassette media used in Horizon® Imagers. This brief pertains to all sizes of DirectVista® Blue-5 Base Film.

Bulletin Details

The various dates are printed on the media supply pack and case box for Horizon media.

Date Format Clarification

Dates printed on all media packaging are in YYYY/MM format. On the case box, this can include the manufacturing date, indicated as  and "Best By" date, indicated as .

"Best By" Date

The media will provide optimal performance when used before the "Best By" date and stored in the conditions indicated below. Optimal media performance may be extended if stored in a cool, dark, dry environment.

Storage for Horizon Media

The media should be stored in the sealed case box. Case boxes should be positioned vertically as indicated by the arrow on the case. Media cases can be stacked two high, vertically.

For high humidity environments, the media should remain

sealed in the original case packaging until just before use. The media should be removed from the case box one pack at a time and the bag inside the case box should be re-sealed around remaining packs.

If the media is cold, be sure to allow the case box to fully reach room temperature before removing individual packs to avoid the risk of moisture condensation on the film.

If the media is taken out of the case box, it should be kept in the cardboard pack box. Media pack boxes should be stored vertically or in the same orientation as in the case box.

Storage Conditions and Requirements

There are three classifications of storage: storage of unprinted media, transportation of unprinted media and archival (extended term) storage of printed media.

- ◆ Unprinted in sealed case box: 0-25° C and 10% - 70% relative humidity
- ◆ Transportation in sealed case box: Maximum 25° C
- ◆ Archival: Extended term storage (20 years) maximum 21° C and 20-30% relative humidity

Technical Support

If problems occur that are not covered by this Technical Brief, contact Codonics Technical Support at any time.

Phone: +1.440.243.1198
Email: support@codonics.com
Website: www.codonics.com

Get it all with just one call
1.800.444.1198

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice.
Patents: www.codonics.com/ip/patents.
Copyright © 2020 Codonics, Inc. Printed in the U.S.A. Part No. 903-417-001.01.



17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
Email info@codonics.com
www.codonics.com

Codonics (Shanghai) Trading Co., Ltd.
317 Xianxia Rd. Building B Unit 1412
Changning District, Shanghai
P.R. China, 200051
86-21-62787701
86-21-62787719 Fax