

# Codonics Grayscale Paper / White Film

Ideal for patient consultations, referring physicians and operating room reference without sacrificing image quality using Codonics exclusive DirectVista® Grayscale Paper



"Paper printing is a boon to radiology, by realizing significant savings in workflow and other soft-costs."

*John C. Litchney, MBA, RT, Administrator for The Cleveland Clinic, Department of Regional Radiology*


"This paper-like media offers many advantages to the end user over traditional film."

*Dr. David Stemerman, Radiologist, Open High-Field MRI and CT of Westchester, NY*

11" x 14"  
28 x 36 cm

14" x 17"  
35 x 43 cm

A and  
A4 Size



- Image quality is far superior to laser and office printers
- Ideal for referring physicians, patient consultations and operating room use
- Letter size allows it to fit conveniently into a referring physician report or patient's chart for quick reference
- No toner, wax or ribbons to ever replace
- No lightbox required
- Provides a great alternative to traditional film

 **CODONICS®**  
We bring the future into focus

# Codonics DirectVista® Grayscale Paper / White Film



DirectVista® Grayscale Paper / White Film is specifically created for medical applications. This high quality media is ideal for referring physicians who prefer to view images in room light without a light box.

Exclusively for Codonics imagers, DirectVista Paper is available in a variety of sizes, including 14"x17" and 11"x 14". Plus, A and A4 size output fits conveniently into a chart so that the image can follow the patient and doctors can easily flip to it, like a note or report.

DirectVista Paper/White Film produces stunning grayscale prints with outstanding contrast ratios. There is no toner, wax or ribbon to replace – ever. It's the ideal media for patient and referral copies as well as archiving.




DirectVista Grayscale Paper/White Film has a heavy, photo-weight base material that resists water and tearing. Since there is no ink or toner, images never smear or transfer onto hands or other prints.

## Specifications:

- Compatible with the Horizon family of multi-media imagers
- Continuous tone grayscale rendering
- Dmin < 0.10 OD, Dmax > 2.00 OD (transmissive)
- Thickness ~ 0.2mm
- Complies with ISO 18911 standards for storage of photographic media
- Fully archivable for a minimum of 10 years with proper storage

## Available Sizes:

- 14" x 17" (35 x 43 cm)
- 11" x 14" (28 x 35 cm)
- A (8.5" x 11")
- A4 (210mm x 297mm)



We bring the future into focus

Call Codonics today at **800.444.1198** or visit [www.codonics.com](http://www.codonics.com) for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents: [www.codonics.com/ip/patents](http://www.codonics.com/ip/patents) Copyright©2006-2022 Codonics Inc. 02/2022



17991 Englewood Drive  
Middleburg Heights, OH 44130 USA  
+1.440.243.1198  
+1.440.243.1334 Fax  
Email [info@codonics.com](mailto:info@codonics.com)  
[www.codonics.com](http://www.codonics.com)